


SUMMARY OF NFSTC® STUDY

VOLUME (ML) OF DECONTAMINANT *Per 20 ML of Blood*

100 ML

200 ML

DNA REDUCTION RESULTS



QTY (NG/ML) OF DNA REMAINING

VOLUME (ML) OF DECONTAMINANT <i>Per 20 ML of Blood</i>	100 ML		200 ML	
	STERIPLEX ^{SD}	0.46	0.81	0.35
BLEACH SOLUTION	5.45	7.42	7.30	1.78
CONTROL <i>NO DECONTAMINANT</i>	12.60	13.50	15.70	12.60